



Giving Physicians More Time to Care System Requirements

Lytec 2010 Practice Management

Workstation system: Equivalent to Core 2 Duo 2.4 processor or higher, 4GB of available hard disk space, 2GB RAM, 32-bit operating system only, and DVD drive

Workstation operating system (required): Windows® XP Professional SP2 or higher, Windows Vista™ Business, 32-bit operating system only

While Lytec now supports Windows 7 and Vista, Revenue Management will only run with Windows XP Professional 32 bit only.

Network: Gigabit full-duplex 10/100/1000 Ethernet router/switch using CAT5e (or CAT6) Ethernet cable. Consult with a licensed network vendor for proper wiring for your business computing environment. Wireless networks are not recommended or supported.

Internet: High-speed cable, DSL or satellite internet service

Server system: Equivalent to: Xeon 2.66Ghz processor or higher, 6GB of available hard disk space, 4GB RAM, and DVD drive

Server operating system (required): Client Server edition: Windows® 2003 or 2008 Server (Enterprise or Standard edition) with latest service packs. Multi-user/ Professional editions: Windows XP Professional SP2 or higher, Windows Vista Business, 32-bit operating system only

Lytec MD EMR

Workstation system: Equivalent to Core 2 Duo 2.4 GHz or higher 2GB of RAM (3GB of RAM for Vista Business) 32GB System Drive or higher 1024x768 Video, 1GB Network Card Windows XP Professional SP2 1500 VA+ UPS Battery Antivirus Solution

Server system: Equivalent to Quad Core 2.66Ghz processor or higher or Xeon 2.66Ghz processor or higher, 4GB RAM, 32GB+ System Drive, 500GB+ Data Drive RAID controller preferred, 1Gbps network card, 1024x768 video, DVD Optical Drive, Tested Backup System, 1500VA+ UPS Battery . Must have separate operating system & data partitions. Data partition must be a non-conflicted network drive and drive letter changed to P:. Domain infrastructure required.

Server operating system
Windows 2003 Server or 2008 Server
32 Bit Only.

Windows Vista Users:
Lytec 2010 currently functions with the Microsoft Windows Vista Business Edition 32 Bit ONLY. Lytec will not function with Vista Ultimate, Home, Home Premium or 64 Bit.

Printer Recommendations:
HP LaserJet 1300 or higher. Microwize will not support any other printers.



Learn more about Lytec at www.microwize.com/lytec